

LONG-NOSE MONOPOLAR FORCEPS - 20 cm



Product Code	30612
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1711053
CND	L1204
EAN/UPC	8023279306125
PARAF	903152902

Description

Insulated Monopolar Instruments:
Bayonet monopolar forceps - 20 cm (7 1/2")

Monopolar forceps provided with a 4 mm connection for standard monopolar cables. It can be used with all GIMA, LED, Erbe and many other electrosurgical units on the market.