


MONOPOLAR SCISSORS - 18 cm



Product Code	30613
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1711053
CND	L1204
EAN/UPC	8023279306132
PARAF	903152799

Description

Insulated Monopolar Instruments:
Monopolar scissors - straight - 18 cm (7")

Monopolar scissors provided with a 4 mm connection for standard monopolar cables. It can be used with all GIMA, LED, Erbe and many other electrosurgical units on the market.