

EU BIPOLAR CABLE FOR MB 80D-120D-160D



Product Code	30617
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No
RDM (NSIS)	2386379
CND	L1899
EAN/UPC	8023279306170

Description

Quality highly flexible silicone bipolar cables 3 m long (10").
Autoclavable at 134°C for minimum 20 times.
Recommended EtO sterilization for max extended life.
Suitable for most generators, with connector for European or US forceps.
Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, PL, CZ, SK, SI, HR, RO, GR, Arabic, JP.

Connector: EU
Compatibility brand/model:
GIMA: MB 80 D, MB 120 D, MB 160 D
LED: SURTRON 80, 120, 160