

US ADSON FORCEPS 12,1 cm



Product Code	30620
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No
RDM (NSIS)	1711035
CND	L1204
EAN/UPC	8023279306200

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: ADSON
- Connector: US
- Length: 12.1 cm (4¾")
- Application: GENERAL/PLASTIC
- Tip size: 1 mm round points