



Product Code	30621
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822650
CND	L1204
EAN/UPC	8023279306217

## Description

### UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: SEMKIN
- Connector: US
- Length: 14 cm (5½")
- Application: GENERAL/PLASTIC
- Tip size: 0.5 mm round points