

# US CUSHING FORCEPS 15.9 cm



Product Code	30622
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1711035
CND	L1204
EAN/UPC	8023279306224

## Description

### UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: CUSHING
- Connector: US
- Length: 15.9 cm (6¼")
- Application: GENERAL/PLASTIC
- Tip size: 1 mm round points