

EU BIPOLAR CABLE FOR MB 240-380



| | |
|---------------|----------------|
| Product Code | 30632 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | Yes |
| CH-REP | No |
| RDM (NSIS) | 2386379 |
| CND | L1899 |
| EAN/UPC | 8023279306323 |

Description

Quality highly flexible silicone bipolar cables 3 m long (10").
Autoclavable at 134°C for minimum 20 times.
Recommended EtO sterilization for max extended life.
Suitable for most generators, with connector for European or US forceps.
Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, PL, CZ, SK, SI, HR, RO, GR, Arabic, JP.

Connector: EU

Compatibility brand/model: GIMA: MB 240, MB 380, VALLEYLAB, CONMED, CODMAN, OLYMPUS