


US BIPOLAR CABLE for MB 120F-200D/F-250-300D-400-400D



| | |
|---------------|----------------|
| Product Code | 30637 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2386379 |
| CND | L1899 |
| EAN/UPC | 8964000208816 |
| GMDN | 10228 |

Description

Quality highly flexible silicone bipolar cables 3 m long (10").
Autoclavable at 134°C for minimum 20 times.
Recommended EtO sterilization for max extended life.
Suitable for most generators, with connector for European or US forceps.
Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, PL, CZ, SK, SI, HR, RO, GR, Arabic, JP.

Connector: US
Compatibility brand/model:
GIMA: MB 120 F, MB 200 F, MB 200 D, MB 200T, MB 300 D, MB 400 D, MB 300T, MB 400T
LED: SURTRON 120, 200, 200T, 300, 300T, 400, 400T
VALLEYLAB, EMC, DOLLEY