



Product Code	30638
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No
RDM (NSIS)	2386379
CND	L1899
EAN/UPC	8964000208786
GMDN	10228

## Description

Quality highly flexible silicone bipolar cables 3 m long (10").  
 Autoclavable at 134°C for minimum 20 times.  
 Recommended EtO sterilization for max extended life.  
 Suitable for most generators, with connector for European or US forceps.  
 Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, PL, CZ, SK, SI, HR, RO, GR, Arabic, JP.

Connector: EU  
 Compatibility brand/model: ERBE, SELECT, DOWNS, SIEMENS