

EU BIPOLAR CABLE for MB 120F-200D/F-250-300D-400-400D



|               |                |
|---------------|----------------|
| Product Code  | 30640          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No  |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2386379       |
| CND        | L1899         |
| EAN/UPC    | 8964000208809 |

Description

Quality highly flexible silicone bipolar cables 3 m long (10").  
Autoclavable at 134°C for minimum 20 times.  
Recommended EtO sterilization for max extended life.  
Suitable for most generators, with connector for European or US forceps.  
Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, PL, CZ, SK, SI, HR, RO, GR, Arabic, JP.

Connector: EU  
Compatibility brand/model:  
GIMA: MB 120 F, MB 200 F, MB 200 D, MB 200T, MB 300 D, MB 400 D, MB 300T, MB 400T  
LED: SURTRON 120, 200, 200T, 300, 300T, 400, 400T  
VALLEYLAB, EMC, DOLLEY