



Product Code	30640
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2386379
CND	L1899
EAN/UPC	8964000208809

## Description

Quality highly flexible silicone bipolar cables 3 m long (10").  
Autoclavable at 134°C for minimum 20 times.  
Recommended EtO sterilization for max extended life.  
Suitable for most generators, with connector for European or US forceps.  
Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, PL, CZ, SK, SI, HR, RO, GR, Arabic, JP.

Connector: EU

Compatibility brand/model:

GIMA: MB 120 F, MB 200 F, MB 200 D, MB 200T, MB 300 D, MB 400 D, MB 300T, MB 400T

LED: SURTRON 120, 200, 200T, 300, 300T, 400, 400T

VALLEYLAB, EMC, DOLLEY