

EU BINNER FORCEPS 18 cm - straight



Product Code	30646
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822897
CND	L1204
EAN/UPC	8964000209899

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: BINNER
- Connector: EU
- Length: 18 cm (7")
- Tip style: STRAIGHT
- Application: ENT
- Tip size: -