

EU BINNER FORCEPS 19 cm - angled 45°



Product Code	30648
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No
RDM (NSIS)	1132406
CND	L1204
EAN/UPC	8023279306484

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: BINNER
- Connector: EUROPEAN
- Length: 19 cm (7 3/8")
- Tip style: ANGLED 45°
- Application: ENT
- Tip size: -