

# EU STRAIGHT FORCEPS 20 cm - angled 1 mm point



Product Code	30653
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822897
CND	L1204
EAN/UPC	8964000450147
GMDN	37020

## Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: STRAIGHT
- Connector: EUROPEAN
- Length: 20 cm (7<sup>3</sup>/<sub>4</sub>"
- Tip style: ANGLED
- Application: GENERAL
- Tip size: 1 mm round points