

EU BAYONET FORCEPS 20 cm - 1 mm point



Product Code	30655
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822897
CND	L1204
EAN/UPC	8964000450192
GMDN	37020

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: BAYONET
- Connector: EUROPEAN
- Length: 20 cm (7³/₄"
- Tip style: STRAIGHT
- Application: GENERAL
- Tip size: 1 mm round points