

EU STRAIGHT FORCEPS 20 cm - 2 mm point



Product Code	30658
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822109
CND	L1204
EAN/UPC	8964000451090
GMDN	37020

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: STRAIGHT
- Connector: EUROPEAN
- Length: 20 cm (7³/₄"
- Tip style: STRAIGHT
- Application: GENERAL
- Tip size: 2 mm round points