

EU STRAIGHT FORCEPS 18 cm - 0.3 mm point



Product Code	30661
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822857
CND	L1204
EAN/UPC	8023279306613
PARAF	903153144
GMDN	37020

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: STRAIGHT
- Connector: EUROPEAN
- Length: 18 cm (7")
- Tip style: STRAIGHT
- Application: GENERAL
- Tip size: 0.3 mm round points