

# EU STRAIGHT FORCEPS 20 cm - 0.3 mm point



Product Code	30662
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1967540
CND	L1204
EAN/UPC	8023279306620
PARAF	903153157

## Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: STRAIGHT
- Connector: EUROPEAN
- Length: 20 cm (7<sup>3</sup>/<sub>4</sub>"
- Tip style: STRAIGHT
- Application: GENERAL
- Tip size: 0.3 mm round points