



Product Code	30664
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822857
CND	L1204
EAN/UPC	8023279306644

## Description

### UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: ADSON
- Connector: EUROPEAN
- Length: 12,1 cm (6")
- Tip style: STRAIGHT
- Application: GENERAL
- Tip size: 1 mm round points