

EU JEWELER FORCEPS 11.5 cm - angled 0.25 mm point



Product Code	30666
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1967540
CND	L1204
EAN/UPC	8023279306668
PARAF	903153106

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: JEWELER 5
- Connector: EUROPEAN
- Length: 11.5 cm (4½")
- Tip style: ANGLED
- Application: OPHTHALMIC/MICRO
- Tip size: 0.25 mm round points