

# EU STRAIGHT FORCEPS 15 cm - angled 0.5 mm point



Product Code	30668
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1967540
CND	L1204
EAN/UPC	8023279306682
PARAF	903153120

## Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: STRAIGHT
- Connector: EUROPEAN
- Length: 15 cm (6")
- Tip style: ANGLED
- Application: OPHTHALMIC/MICRO
- Tip size: 0.5 mm round points