

EU MC-PHERSON FORCEPS 10,2 cm - straight 0.5 mm point



Product Code	30669
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1967540
CND	L1204
EAN/UPC	8023279306699

Description

UNIVERSAL HIGH QUALITY INSULATED BIPOLAR FORCEPS

- Model: Mc PHERSON
- Connector: EUROPEAN
- Length: 10.2 cm (4")
- Tip style: STRAIGHT
- Application: OPHTHALMIC/MICRO
- Tip size: 0.5 mm round points