

# SPATULA LAPAROSCOPIC ELECTRODE - 47 cm



Product Code	30685
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CND

K02030101

## Description

### LAPAROSCOPIC STERILE ELECTRODES - LONGER HANDLE

Spatula electrode - 47 cm

Long handle gives the surgeon a better grip and makes the instrument more balanced.

Laparoscopic electrodes without suction can be used with any electrosurgery unit connecting cable 4 mm (30609) to monopolar output.