

LOOP ELECTRODE 15 x 8 mm - autoclavable



Product Code	30687
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1710953
CND	K02030101
EAN/UPC	8023279306873
PARAF	903152015

Description

LOOP electrode

- Ø 2.4 mm
- autoclavable at 134°C
- non sterile
- Shaft length: 13 cm
- Electrode width: 15 mm
- Electrode depth: 8 mm