

ETHICON ETHILON MONOFILAMENT SUTURES - gauge 5/0 needle 13 mm



Product Code	22307
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2561758
CND	H010201010101
EAN/UPC	8023279223071
GMDN	38000

Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 698H
- Gauge: 5/0
- Metric: 1.0
- Colour: black
- Circle: 3/8
- Needle lenght: 13 mm
- Needle: P-3
- Suture lenght: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.