



Product Code	30700
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2783315
CND	Z12010901
EAN/UPC	8033003260296
GMDN	35029

Description

ELECTROCAUTERY - 1,200°C - fine tip - box of 10

Electrocautery is a sterile bipolar battery operated disposable device for the cauterization of small vessels. With a fine tip required in most of the surgeries.

Long life internal battery (60 min.)

Length: 18 cm

Made in Italy.

Applications:

- Cardiovascular, thoracic surgery
- Proctology
- Gynaecology
- Neurosurgery
- Orthopedic
- Urology
- Obstetrics
- Otolaryngology
- Plastic and reconstructive
- Ophthalmology
- Dermatology
- Dental surgery

Technical Specifications

- HTC (High Temperature Cautery)
- Tip shape: Fine
- Tip length: 28mm
- Temperature °C/°F: 1200/2200
- Weight: 65g