

# ELECTROCAUTERY 1200°C - large tip



Product Code	30701
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2783310
CND	Z12010901
EAN/UPC	8033003262979
GMDN	35029

## Description

ELECTROCAUTERY - 1,200°C - large tip - box of 10

Electrocautery is a sterile bipolar battery operated disposable device for the cauterization of small vessels.

Long life internal battery (60 min.)

Length: 18 cm

Made in Italy.

Applications:

- Proctology
- Neurosurgery
- Orthopedic
- Otolaryngology
- Plastic and reconstructive

## Technical Specifications

- HTC (High Temperature Cautery)
- Tip shape: Large
- Tip length: 28mm
- Temperature °C/°F: 1200/2200
- Weight: 65g