

ELECTROCAUTERY 1200°C - long shaft - fine tip



Product Code	30702
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2783275
CND	Z12010901
EAN/UPC	8033003260104
GMDN	35029

Description

LONG SHAFT ELECTRO CAUTERY - 1,200°C - fine tip - box of 10

Electrocautery is a sterile bipolar battery operated disposable device for the cauterization of small vessels. With the tip fixed on a 12 cm long shaft.

Precise cauterization even in difficult access sites.

Long life internal battery (60 min.)

Made in Italy.

Applications:

- Proctology
- General surgery
- Surgical gynaecology
- Neurosurgery
- Otolaryngology

Technical Specifications

- HTC (High Temperature Cautery)
- Tip shape: Fine
- Tip length: 125mm
- Temperature °C/°F: 1200/2200
- Weight: 78g