

# EMERGENCY ELECTROCAUTERY 800°C - thick tip



Product Code	30703
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2783305
CND	Z12010901
EAN/UPC	8033003260173
GMDN	35029

## Description

EMERGENCY ELECTROCAUTERY - 800°C - thick tip - box of 10

Electrocautery is a sterile bipolar battery operated disposable device for the cauterization of small vessels. With "thick tip" for the perforation of nails in case of haematoma.

Long life internal battery (60 min.)

Made in Italy.

Applications:  
• Emergency ward

## Technical Specifications

- HTC (High Temperature Cautery)
- Tip shape: Thick
- Tip length: 28mm
- Temperature °C/°F: 800/1200
- Weight: 65g