



Product Code	30704
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2783277
CND	Z12010901
EAN/UPC	8033003260098
GMDN	35029

Description

ELECTROCAUTERY - 600°C - fine tip - box of 10

Electrocautery is a sterile bipolar battery operated disposable device for the cauterization of small vessels. With lower temperature 600°C. Especially designed for ophthalmology.

Long life internal battery (60 min.)

Made in Italy.

Applications:
• Ophthalmology

Technical Specifications

- LTC (Low. Temp. Cautery)
- Tip shape: Fine
- Tip length: 28mm
- Temperature °C/°F: 600/1100
- Weight: 45g