

GIMANORD LED PLUS MAGNIFYING LIGHT tunable white - trolley



Product Code	30715
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2950258
CND	Z12010702
EAN/UPC	8023279307153

Description

GIMANORD PLUS TUNABLE WHITE - on trolley

They feature a biconvex magnifying lens in optical glass with Ø 120 mm and 3 dioptries. The cold LED light allows prolonged use, without causing any eyesight problem. Polycarbonate screen for protection of the light source. Slim and straightforward reflector. Lamps are lightweight and handling. Suitable for dermatology use. Medical device Class I.

Touch control button

It allows to adjust light intensity from 4% to 100% and colour temperature from 2,700 to 6,000 °K with a smooth touch.

Interchangeable plugs

The transformer, integrated in the plug for a lighter and smooth movements of the arm structure, is supplied with interchangeable plugs (UK, US, EU).

QR Manual: GB, FR, IT, ES, PT

Made in Italy.

Technical Specifications

Light intensity at 0.5 m (Ec) 2,000 lux
Colour temperature 2,700-6,000 K
Colour rendering index (CRI) 95 Ra
R9 =90
Operating voltage 100-240 V - 50/60 Hz
Absorbed power 16 W
Light head diameter 23 cm
Average life >60.000 hours
Norms IEC-60601-1, IEC-60601-2-41, EU 2017/745