



Product Code	30716
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2391800
CND	Z12010702
EAN/UPC	8023279307160

Description

GIMANORD PLUS - with table clamp

This model features a biconvex magnifying lens in optical glass with Ø 120 mm. The cold LED light allows prolonged use, without causing any eyesight problem. Polycarbonate screen for protection of the light source. Suitable for dermatology use. Medical device Class I.

Design

The new design makes the lamp extremely easy to use and manoeuvre, with an additional rotation of the reflector.

Touch control

Touch button allows to adjust light intensity from 4% to 100% with a smooth touch.

Interchangeable plugs

The transformer, integrated in the plug for a lighter and smooth movements of the arm structure, is supplied with interchangeable plugs (UK, US, EU).

QR Manual: GB, FR, IT, ES, PT

Made in Italy.

Technical Specifications

Light intensity at 0.5 m (Ec) 2,250 lux
 Colour temperature 5,370 K
 Colour rendering index (CRI) 95 Ra
 R9 =90
 Operating voltage 100-240 V - 50/60 Hz
 Absorbed power 16.5 W
 Light head diameter 23 cm
 Average life >60.000 hours
 Norms IEC-60601-1, IEC-60601-2-41, EU 2017/745