



Product Code	30736
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279307368
---------	---------------

PARAF	903113746
-------	-----------

Description

DERMALED MAGNIFYING LAMP

Cold light LED lamp with 3 dioptries magnifying glass lens and brightness adjustment. It comes with a high stability 5 wheels base and an articulated arm for easy movement and fixing of the lamp.

Manual in GB, IT.

Technical Specifications

ype of light: LED

Voltage: 7 W

Dioptries: 3 (large glass)

Colour temperature: 7,500-9,500 °K

Light intensity at 20 cm: 2,400 lux

Operating voltage: 100-240 V - 50/60 Hz

Net/Gross weight: 12.8-15 kg