



Product Code	30737
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1816543
CND	Z12010701
EAN/UPC	8436562860073

Description

LUPA LED MAGNIFYING LIGHT - trolley

Smartly designed white cold light LED lamp with 5-dioptres magnification. Generates a light without shadow. 5 wheels base and articulated arm allow easy movement.

Manual in GB, FR, IT, ES, PT, DE.

Made in Spain.

Height 90 cm (arm) + 95 cm (trolley).
Weight: 11.8 Kg (8.8 Kg trolley + 3 Kg arm).

Technical Specifications

Type of light: 60 LED - HF 6 W
 Dioptres: 5
 Lens: Ø 12.7 cm
 Colour temperature: 6,500 °K
 Operating voltage: 110-230 V - 50/60 Hz
 Class I
 Marks: 93/42/CEE
 Lamp life: 40,000 hours