



Product Code	30739
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8436562860073

## Description

LUPA LED MAGNIFYING LIGHT - desk with clamp

Smartly designed white cold light LED lamp with 5-dioptres magnification. Generates a light without shadow.

Manual in GB, FR, IT, ES, PT, DE.

Made in Spain.

## Technical Specifications

Type of light: 60 LED - HF 6 W

Dioptres: 5

Lens: Ø 12.7 cm

Colour temperature: 6,500 °K

Operating voltage: 110-230 V - 50/60 Hz

Class I

Marks: 93/42/CEE

Lamp life: 40,000 hours