

LUXIFLEX HALOGEN LIGHT 50,000 Lux - trolley



Product Code	30754
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1984336
CND	Z12010702
EAN/UPC	8436562860004

Description

LUXIFLEX on trolley - halogen

Warm light luminaire, versatile with compact design.
Designed for examination procedures that require great performance. Easy use, thanks to its compact head withØ 8 cm and its switch located at thebase of the lamp.
Class II.

Manual in: GB, FR, IT, DE, ES, PT

Height 70 cm (arm) + 92 cm (trolley)

Made in Spain.

Technical Specifications

Type of lighting: halogen
Consumption (Watt): 35
Life of the bulbs (hours): > 5,000
Light intensity at 50 (lux): > 35,000
Light intensity at 1 m (lux): 12,000
Colour temperature at 5% at 1 m (K): 3,000
Colour rendering index (Ra)%: 100
Light field size at 50 cm (Ø mm): 150
External Ø of the reflector (mm): 80
Colour: white
Norms: EN 60601-1-2
Marks: EU 2017/745
Weight (kg): -
Operating voltage V: 230
Operating Hz: 50/60