



Product Code	30754
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I
UK-REP	No
CH-REP	No
RDM (NSIS)	1984336
CND	Z12010702
EAN/UPC	8436562860004

Description

LUXIFLEX on trolley - halogen

Warm light luminaire, versatile with compact design.

Designed for examination procedures that require great performance. Easy use, thanks to its compact head with \varnothing 8 cm and its switch located at the base of the lamp.

Class II.

Manual in: GB, FR, IT, DE, ES, PT

Height 70 cm (arm) + 92 cm (trolley)

Made in Spain.

Technical Specifications

Type of lighting: halogen
 Consumption (Watt): 35
 Life of the bulbs (hours): > 5,000
 Light intensity at 50 (lux): > 35,000
 Light intensity at 1 m (lux): 12,000
 Colour temperature at 5% at 1 m (K): 3,000
 Colour rendering index (Ra)%: 100
 Light field size at 50 cm (\varnothing mm): 150
 External \varnothing of the reflector (mm): 80
 Colour: white
 Norms: EN 60601-1-2
 Marks: EU 2017/745
 Weight (kg): -
 Operating voltage V: 230
 Operating Hz: 50/60