

PORTABLE LED LIGHT SOURCE 200.000 LUX



Product Code	30796
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2567717
CND	Z12020402
EAN/UPC	6975465800066

Description

PORTABLE LED LIGHT SOURCE - 200,000 LUX

Powerful portable LED light source mainly used in endoscopic examinations, such as ENT. Its button switch can adjust the light on 3 levels. 5W LED with condenser system.

Supplied with 2 adaptors for Storz and Wolf endoscopes and 1 adaptor for fiber optic cable, lithium battery, charger and cable. GB, IT manual.

Needs 30832 for use with GIMA F.O. cables

Size: 136x24x24 mm

Technical Specifications

Illumination: > 200,000 lux
Colour temperature: 5,700/6,500 k
Bulb power: 5 W LED
Rechargeable battery: li-ion, 2,600mA
Weight: 90 g
Power supply: AC 100-240 V
Charging time: 3 hours
CRI: 75-92