



Product Code	30803
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2215903
CND	Z129004
EAN/UPC	8023279308037

Description

GIMALED 200 K LIGHT SOURCE - 16 LED

Compact, lightweight, powerful light source for endoscopic use.
A low noise fan provides the necessary cooling inside the metallic structure of the light source.
Touch switch enables a variable light intensity.

Weight: 3 kg
Optional connectors allow use with STORZ, WOLF, OLYMPUS, PENTAX, FUJINON cables

User manual in GB, ES, FR, IT, PT.

Made in Italy.

Technical Specifications

Input for: GIMA, SCHOLLY and LUT cables
Illumination at 2.5 cm: 147,000 Lux \pm 10 %
Colour temperature: 5000 K \pm 5%
CRI (Colour Rendering Index): 89 Ra
Max power absorption: 34 VA
Power supply: 100-240V, 50/60 Hz
LED life: 60,000 hours
Dimension: 27x17x h 11 cm
Norms: IEC 60601-1, IEC 60601-2-18 Class II