

# F.O.CABLE 4.8x1800 mm - without adapters



Product Code	30818
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Not required
CH-REP	Not required
RDM (NSIS)	2799312
CND	V9099
EAN/UPC	8023279308181
PARAF	903145237

## Description

High transmission surgical light cable, armoured but highly flexible with unique features:

- no burn outs
- strain relieved and protected by silicone hose
- steam autoclavable with warranty (12 months)
- light source end temperature resistant up to 134°C
- kink resistant
- no electrical interferences with endoscopic surgery units
- surgical stainless steel end fittings resist wear and corrosion
- stainless steel adaptors enables our cable to be used with other manufacturers halogen light sources/instruments