



Product Code	30848
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279308488
---------	---------------

Description

SMART MAGNIFYING GLASSES 2.5X

Spectacle frame with magnifying glasses allows the physician to have a magnified view of the visual field.
Lightweight, it folds up and can be worn over spectacles.
2.5X zoom lens can be tilted to achieve optimum focus.