

RIESTER RI-FOCUS 6091 LED HEADLIGHT



Product Code	30881
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1961018
CND	Z129004
EAN/UPC	4045396135680

Description

High performance LED headlight with battery compartment for fast and easy change of battery. Flexible hinge for a comfortable one-hand use of the LED illumination body.

- more efficient diagnosis due to white and high performance LED with up to 10,000 lux and a 50,000 hours lifetime
- wireless headband with built-in battery compartment
- particularly comfortable with balanced headband and one-hand operation
- focusable LED lamp at a distance of 40 cm approx. Ø 40 to 200 mm infinitely focusable
- white LED with 6000°K, 140 lumen
- significantly less heating and energy loss
- rechargeable battery operation: approx. 90 min with fully charged batteries
- rechargeable AAA batteries and plug in 230V charger included

On/off switch on the battery compartment
Removable and washable inner padding
Weight: 210 g without batteries

Made in Germany.

Technical Specifications

- Distance: 400 mm/200 mm
- Lux max: 10,000/23,000
- Circle of light max Ø: 200 mm/120 mm
- Circle of light min Ø: 40 mm/30 mm