

HEINE ML4 LED HQ HEAD LIGHT PLUG IN TRANSFORMER J.008.31.413



Product Code	30884
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2750133
CND	Z129004
EAN/UPC	4053755180318

Description

HEINE ML4 LEDHQ HEADLIGHT PLUG IN TRANSFORMER- J.008.31.413

- LED in HQ; defining optimal light intensity, brightness and homogeneity from edge to edge, and colour rendering (4,500 k, CRI >90, R9 >75). This ensures true colour tissue for the most accurate diagnosis.
- Coaxial illumination: ensures a completely shadow-free image and allows for excellent illumination of difficult to see areas
- Lithium Battery: for up to 8.5 hours of continuous working
- LED Thermal Management: guarantees consistent quality and performance throughout the LED's lifetime (50,000 hours)
- Adjustable illumination spot size: 30 to 80 mm spot size range (at 420 mm working distance) to fit any examination situation
- Comfortable and secure fit: multiple adjustment points and soft padding for comfort and stability even during long procedures
- Illumination Intensity Control: for precise brightness setting

Made in Germany.

Technical Specifications

Power supply: mPack® Rechargeable Battery Pack with Li-ion battery and direct mains transformer
Illumination: 90,000 lux at a working distance of 180 mm
LED operating life: 50,000 hours
Spot adjustment: 30 mm to 80 mm illuminated field diameter (at 420 mm)
Colour temperature: 5,500°K