

ETHICON PERMA-HAND SILK SUTURES - gauge 4/0 needle 17 mm RB-1 - braided



Product Code	22309
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2471107
CND	H010202020301
GTIN	10705031044562

Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: K871H
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: 1/2
- Needle length: 17 mm
- Needle: RB-1
- Suture length: 75 cm
- Minimum order: box of 36 pieces

Perma-Hand silk suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures.