

STYLE BINOCULAR LOUPE 3.5x 340mm



Product Code	30914
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2624944
CND	Z12119099
EAN/UPC	8023279309140
GTIN	08023279309140
GMDN	35843

Description

STYLE BINOCULAR LOUPE

Line of binocular loupes with same features of traditional line but mounted on a extra light, fashionable frame.

- Style binocular Loupe
- Magnification: 3.5X
- Working distance: 340 mm
- Field of view Ø: 50 mm

Lightweight, ergonomically superior Galilean-style flip-up loupes.

Can be used independently or mounted on 30890, 30940.

Multi-coated and high grade glass lenses provide clean and sharp image.

Eyepiece with wide field of view.

Deep focus to see object stereoscopically.

Loupes angles can be easily adjusted between 0° and 40° to minimize fatigue even after prolonged use.

6 mm up and down displacement.

Adjustable pupil distance to fit different face types.

Weight: 111 g

Multilingual manual: GB, FR, IT, ES, PT, GR.

All binocular loupes are supplied in a case