



Product Code	30932
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	96158
CND	Z12019080
EAN/UPC	4053755118274
UMDNS	15-652
GMDN	32692

Description

HEINE® S-FRAME® WITH BINOCULAR LOUPES

Heine binocular loupes

- Magnification: 4x
- Working distance: 340 mm
- HEINE code: C-000.32.431

- achromatic optics: crisp, bright image with colour-correction.
- multi-coated, scratch-resistant silica glass lenses
- produce crisp, clear images.

High magnification with a large field of view and excellent depth of field.

- Water resistant
- Flexibility for precise adjustment: independently adjustable Pupillary Distance controls ensure precise alignment of the optics, for a stereoscopy view.
- i-View loupe mount; provides any angle of view and the possibility flips up and down the optics to exclude or reinclude the magnification during any procedure.
- Different magnification and working distance

Standard accessories

Set content:

- optics
- S-frame with i-View loupe mount
- detachable protective lenses
- 2 sterilisable swivel levers
- retaining cord
- carrying case