



Product Code	30936
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2269745
CND	Z129004
EAN/UPC	8023279309362

Description

GIMA 50 KL SURGICAL LED HEADLIGHT

Latest generation of portable headlight system offering high quality performance.

- spot size adjustable by inner mini light grating
- precise light imaging by multi optical lens
- integrated cordless design

Supplied with battery and battery charger.

Multilingual manual: GB, FR, IT, ES, PT, DE, DK, SE, PL, HU, SK, RO, HR, GR.

Technical Specifications

Illumination source: 5W LED
 Illumination intensity: 10-50,000 lux
 Colour temperature: 6,000K
 Spot size: 50-100 at 20-50 cm
 Lamp working life: 50,000 hours
 Battery: 3.7V Lithium Battery
 Battery run time: 2-12 hours
 Battery charging time: 2 hours
 Input voltage: 100-240VAC
 Weight: 200 g