

# RIESTER SUPERVU GALILEAN BINOCULAR LOUPE

## 2.0x 340mm - 3960-620



Product Code	30950
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 4045396160606

### Description

#### SUPERVU GALILEAN BINOCULAR LOUPES

Riester code: 3960-620  
Magnif.: 2.0 X  
Working distance: 340 mm  
Field of view Ø: 114 mm

They offer 2.0 X magnification which allows greater distance between the user and the patient.  
They provide a carefully calculated blend of comfort, clarity and field of view.

- excellent optical features with respect to the field view ( $\pm 10\%$ )
- movement whilst maintaining focal clarity is possible due to superior depth of field.

The lightweight aluminium bar with individual pupillary adjustment provides the mounted loupes with stability and optimal focal alignment, eliminating headache and eyestrain

- ideal for ocular and dental examinations as well as plastic surgery and veterinary professionals.

99 g (including frame).