

RIESTER SUPERVU HI-RES GALILEAN BINOCULAR LOUPE 3.0x 340mm - 3970-630



Product Code	30960
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2572657
CND	Y210321
EAN/UPC	4045396160705

Description

SUPERVU HIRES GALILEAN BINOCULAR LOUPES

Riester code: 3970-630
Magnif.: 3.0 X
Working distance: 340 mm
Field of view Ø: 68 mm

They provide a high-resolution loupe, which far surpasses any previous standards.

This 5-element loupe offers triple magnification with a 15% increased field of view, a bright image and a visually flat field from edge to edge.

- excellent optical features with respect to the field view ($\pm 10\%$)
- 3.0 X magnification
- excellent distinction, stunning performance as well as great looks
- ideal for most fields of surgery, surgical assistants as well as ophthalmologists