

RIESTER SUPERVU HI-RES GALILEAN BINOCULAR LOUPE 3.0x 420mm - 3971-630



| | |
|---------------|----------------|
| Product Code | 30961 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2572657 |
| CND | Y210321 |
| EAN/UPC | 4045396160712 |

Description

SUPERVU HIRES GALILEAN BINOCULAR LOUPES

Riester code: 3971-630
Magnif.: 3.0 X
Working distance: 420 mm
Field of view Ø: 82 mm

They provide a high-resolution loupe, which far surpasses any previous standards.

This 5-element loupe offers triple magnification with a 15% increased field of view, a bright image and a visually flat field from edge to edge.

- excellent optical features with respect to the field view ($\pm 10\%$)
- 3.0 X magnification
- excellent distinction, stunning performance as well as great looks
- ideal for most fields of surgery, surgical assistants as well as ophthalmologists